



Provox® LaryTube™


Provox LaryTube



Provox LaryTube, Fenestrated



Provox LaryTube with Ring



Product description:

The Provox LaryTube is a tube made of medical grade silicone rubber. The purpose of the device is to create a comfortable and airtight fit between the Provox LaryTube and the tracheostoma, and to provide attachment for devices from the Provox HME System.

The Standard model and the Ring version can be fenestrated so that air can go through the voice prosthesis for voice prosthesis users.

Standard versions – made for use with or without a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClip.

Fenestrated versions – for voice prosthesis users. Can be attached with a Provox TubeHolder or Provox LaryClip.

Ring versions – made for use with or without a voice prosthesis.

Document ID: PF011-01-TechInfo

Edition: 2.0

Manufacturer: Atos Medical AB
Kraftgatan 8
SE-242 35 Hörby, Sweden

Classification: (EU) IIb (Rule 5)
93/42/EEC

Intended Use: The Provox LaryTube is a holder for devices in the Provox HME System intended for vocal and pulmonary rehabilitation after total laryngectomy. For patients with a shrinking tracheostoma it is also used to maintain the tracheostoma for breathing. The Provox LaryTube is intended for single patient use.

Use specifications: **Intended medical indication**
Product for rehabilitation for patients breathing through a tracheostoma.

Intended patient population

Patients of any age.
Cognitive ability, by a clinician judged as sufficient.
Manual dexterity, by a clinician judged as sufficient.

Intended usage

Single patient multiple use, Prescription only.

Intended part of the body/type of tissue applied to or interacted with

Tracheostoma.

Intended user profile

The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.

Intended conditions of use

Environment: Home use (normal daily environment without any or environmental restrictions regarding temperature, moisture etc.).
Outpatient clinic use. Hospital use.
Frequency of use: Continuous use.
Replacement rate: Max use of 6 months. Replacement is performed by the patient, clinician or caregiver.

Contraindications: Provox LaryTube is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheal or tracheostoma tissue.

CE Mark: Yes. Devices are CE-marked.

GMDN code: 12292 (Laryngectomy tube)

Sterilization: Non-sterile

Raw material: LaryTube: Silicone
Ring: Silicone with blue masterbatch

Latex information: Not manufactured with natural rubber latex.

Product Information

| | |
|-------------------------------------|--|
| Biological origin: | The device is not manufactured with materials derived from human or animal source. |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2 °C - 42 °C. |
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. |
| Hazardous components: | None. |
| Expiration date: | 3 years after manufacturing. |
| Packaging: | <p>Provox LaryTube (standard) is packed in a plastic bag of polyethylene. 5 pcs of Provox XtraFlow HME are packed in a plastic bag of polyethylene. The products and instructions for use for Provox LaryTube and Provox XtraHME are packed in a cardboard box.</p> <p>Provox LaryTube (fenestrated) is packed in a plastic bag of polyethylene. 5 pcs of Provox XtraFlow HME are packed in a plastic bag of polyethylene. 1 Provox Brush is packed in a plastic bag of polyethylene. The products and instructions for use for Provox LaryTube and Provox XtraHME and Provox Brush are packed in a cardboard box.</p> <p>Provox LaryTube (with ring) is packed in a plastic bag of polyethylene. 5 pcs of Provox XtraFlow HME are packed in a plastic bag of polyethylene. The products and instructions for use for Provox LaryTube and Provox XtraHME are packed in a cardboard box.</p> |

Devices under Basic UDI-DI: 7331791-LTU-0-000-0002-3E

| REF | Name | UDI-DI |
|------|------------------------------------|----------------|
| 7601 | Provox LaryTube 8/27 | 07331791002076 |
| 7602 | Provox LaryTube 8/36 | 07331791002090 |
| 7603 | Provox LaryTube 8/55 | 07331791002113 |
| 7605 | Provox LaryTube 9/27 | 07331791002137 |
| 7606 | Provox LaryTube 9/36 | 07331791002151 |
| 7607 | Provox LaryTube 9/55 | 07331791002175 |
| 7609 | Provox LaryTube 10/27 | 07331791002199 |
| 7610 | Provox LaryTube 10/36 | 07331791002212 |
| 7611 | Provox LaryTube 10/55 | 07331791002236 |
| 7613 | Provox LaryTube 12/27 | 07331791002250 |
| 7614 | Provox LaryTube 12/36 | 07331791002274 |
| 7615 | Provox LaryTube 12/55 | 07331791002298 |
| 7624 | Provox LaryTube 8/36 with Ring | 07331791002311 |
| 7625 | Provox LaryTube 8/55 with Ring | 07331791002335 |
| 7626 | Provox LaryTube 9/36 with Ring | 07331791002359 |
| 7627 | Provox LaryTube 9/55 with Ring | 07331791002373 |
| 7628 | Provox LaryTube 10/36 with Ring | 07331791002397 |
| 7629 | Provox LaryTube 10/55 with Ring | 07331791002410 |
| 7630 | Provox LaryTube 12/36 with Ring | 07331791002434 |
| 7631 | Provox LaryTube 12/55 with Ring | 07331791002458 |
| 7637 | Provox LaryTube 8/36, Fenestrated | 07331791002472 |
| 7638 | Provox LaryTube 8/55, Fenestrated | 07331791002496 |
| 7640 | Provox LaryTube 9/36, Fenestrated | 07331791002519 |
| 7641 | Provox LaryTube 9/55, Fenestrated | 07331791002533 |
| 7643 | Provox LaryTube 10/36, Fenestrated | 07331791002557 |
| 7644 | Provox LaryTube 10/55, Fenestrated | 07331791002571 |
| 7646 | Provox LaryTube 12/36, Fenestrated | 07331791002595 |
| 7647 | Provox LaryTube 12/55, Fenestrated | 07331791002618 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|-------------------------------|---------------------------|
| Provox Adhesive | 7331791-ADH-0-000-0000-CQ |
| Provox BasePlate Adaptor | 7331791-HME-A-000-0003-F5 |
| Provox Brush | 7331791-VPS-A-000-0003-RR |
| Provox FreeHands HME Cassette | 7331791-HME-0-000-0003-XJ |
| Provox LaryClip | 7331791-LTU-A-000-0001-JT |
| Provox Micron HME | 7331791-HME-0-000-0002-XF |
| Provox ShowerAid | 7331791-ADH-A-000-0000-U8 |
| Provox Swab | 7331791-GEN-A-000-0002-EC |
| Provox TubeBrush | 7331791-GEN-A-000-0001-E9 |
| Provox TubeHolder | 7331791-GEN-A-000-0000-E6 |
| Provox XtraHME | 7331791-HME-0-000-0000-X9 |


Document Approvals
Approved Date: 2023-10-24

| | |
|---|---|
| Task: Approval Task Verdict: Approve | ABDALM Abdallah Almashharawi, Sustaining Engineer (abdallah.almashharawi-atosmedical@coloplast.com) Issuer 23-Oct-2023 10:05:09 GMT+0000 |
|---|---|


| | |
|--|---|
| Task: Final Approval Verdict: Approve | ELIAND Elin Andersson, Associate Design Control & Usability Specialist (elin.andersson-atosmedical@coloplast.com) Technical / Specialist 24-Oct-2023 09:29:44 GMT+0000 |
|--|---|

Provox® LaryTube™


Provox LaryTube



Provox LaryTube, Fenestrated



Provox LaryTube with Ring



Product description:

The Provox LaryTube is a tube made of medical grade silicone rubber. The purpose of the device is to create a comfortable and airtight fit between the Provox LaryTube and the tracheostoma, and to provide attachment for devices from the Provox HME System.

The Standard model and the Ring version can be fenestrated so that air can go through the voice prosthesis for voice prosthesis users.

Standard versions – made for use with or without a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClip.

Fenestrated versions – for voice prosthesis users. Can be attached with a Provox TubeHolder or Provox LaryClip.

Ring versions – made for use with or without a voice prosthesis.

Document ID: PF011-02-TechInfo **Edition:** 2.0

Manufacturer: Atos Medical AB
Kraftgatan 8
SE-242 35 Hörby, Sweden

Classification: (EU) IIb (Rule 5)
93/42/EEC

Intended Use: The Provox LaryTube is a holder for devices in the Provox HME System intended for vocal and pulmonary rehabilitation after total laryngectomy. For patients with a shrinking tracheostoma it is also used to maintain the tracheostoma for breathing. The Provox LaryTube is intended for single patient use.

Use specifications: **Intended medical indication**
Product for rehabilitation for patients breathing through a tracheostoma.

Intended patient population

Patients of any age.
Cognitive ability, by a clinician judged as sufficient.
Manual dexterity, by a clinician judged as sufficient.

Intended usage

Single patient multiple use, Prescription only.

Intended part of the body/type of tissue applied to or interacted with

Tracheostoma.

Intended user profile

The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.

Intended conditions of use

Environment: Home use (normal daily environment without any or environmental restrictions regarding temperature, moisture etc.).
Outpatient clinic use. Hospital use.
Frequency of use: Continuous use.
Replacement rate: Max use of 6 months. Replacement is performed by the patient, clinician or caregiver.

Contraindications: Provox LaryTube is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheal or tracheostoma tissue.

CE Mark: Yes. Devices are CE-marked.

GMDN code: 12292 (Laryngectomy tube)

Sterilization: Non-sterile

Raw material: LaryTube: Silicone
Ring: Silicone with blue masterbatch

Latex information: Not manufactured with natural rubber latex.

Product Information

| | |
|-------------------------------------|--|
| Biological origin: | The device is not manufactured with materials derived from human or animal source. |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2 °C - 42 °C. |
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. |
| Hazardous components: | None. |
| Expiration date: | 5 years after manufacturing. |
| Packaging: | Provox LaryTube is packed in a plastic bag of polyethylene. The product and instructions for use are packed in a cardboard box. |

Devices under Basic UDI-DI: 7331791-LTU-0-000-0002-3E

| REF | Name | UDI-DI |
|--------|------------------------------------|---------------|
| 7601FR | Provox LaryTube 8/27 | 7331791002083 |
| 7602FR | Provox LaryTube 8/36 | 7331791002106 |
| 7603FR | Provox LaryTube 8/55 | 7331791002120 |
| 7605FR | Provox LaryTube 9/27 | 7331791002144 |
| 7606FR | Provox LaryTube 9/36 | 7331791002168 |
| 7607FR | Provox LaryTube 9/55 | 7331791002182 |
| 7609FR | Provox LaryTube 10/27 | 7331791002205 |
| 7610FR | Provox LaryTube 10/36 | 7331791002229 |
| 7611FR | Provox LaryTube 10/55 | 7331791002243 |
| 7613FR | Provox LaryTube 12/27 | 7331791002267 |
| 7614FR | Provox LaryTube 12/36 | 7331791002281 |
| 7615FR | Provox LaryTube 12/55 | 7331791002304 |
| 7624FR | Provox LaryTube 8/36 with Ring | 7331791002328 |
| 7625FR | Provox LaryTube 8/55 with Ring | 7331791002342 |
| 7626FR | Provox LaryTube 9/36 with Ring | 7331791002366 |
| 7627FR | Provox LaryTube 9/55 with Ring | 7331791002380 |
| 7628FR | Provox LaryTube 10/36 with Ring | 7331791002403 |
| 7629FR | Provox LaryTube 10/55 with Ring | 7331791002427 |
| 7630FR | Provox LaryTube 12/36 with Ring | 7331791002441 |
| 7631FR | Provox LaryTube 12/55 with Ring | 7331791002465 |
| 7637FR | Provox LaryTube 8/36, Fenestrated | 7331791002489 |
| 7638FR | Provox LaryTube 8/55, Fenestrated | 7331791002502 |
| 7640FR | Provox LaryTube 9/36, Fenestrated | 7331791002526 |
| 7641FR | Provox LaryTube 9/55, Fenestrated | 7331791002540 |
| 7643FR | Provox LaryTube 10/36, Fenestrated | 7331791002564 |
| 7644FR | Provox LaryTube 10/55, Fenestrated | 7331791002588 |
| 7646FR | Provox LaryTube 12/36, Fenestrated | 7331791002601 |
| 7647FR | Provox LaryTube 12/55, Fenestrated | 7331791002625 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|-------------------------------|---------------------------|
| Provox Adhesive | 7331791-ADH-0-000-0000-CQ |
| Provox BasePlate Adaptor | 7331791-HME-A-000-0003-F5 |
| Provox Brush | 7331791-VPS-A-000-0003-RR |
| Provox FreeHands HME Cassette | 7331791-HME-0-000-0003-XJ |
| Provox LaryClip | 7331791-LTU-A-000-0001-JT |
| Provox Micron HME | 7331791-HME-0-000-0002-XF |
| Provox ShowerAid | 7331791-ADH-A-000-0000-U8 |
| Provox Swab | 7331791-GEN-A-000-0002-EC |
| Provox TubeBrush | 7331791-GEN-A-000-0001-E9 |
| Provox TubeHolder | 7331791-GEN-A-000-0000-E6 |
| Provox XtraHME | 7331791-HME-0-000-0000-X9 |

Document Approvals
Approved Date: 2023-10-24


| | |
|---|---|
| Task: Approval Task Verdict: Approve | ABDALM Abdallah Almashharawi, Sustaining Engineer (abdallah.almashharawi-atosmedical@coloplast.com) Issuer 23-Oct-2023 10:05:08 GMT+0000 |
|---|---|

| | |
|--|---|
| Task: Final Approval Verdict: Approve | ELIAND Elin Andersson, Associate Design Control & Usability Specialist (elin.andersson-atosmedical@coloplast.com) Technical / Specialist 24-Oct-2023 09:27:50 GMT+0000 |
|--|---|

| Justification: | Function: | Electronic signature justification: | Signed: Date (yyyy-mm-dd) - Time (hh:mm): |
|----------------|-----------|-------------------------------------|---|
| Issued: | QA | Carolina Johansson - SEHRBJNC | 2022-04-12 - 14:46 |
| Reviewed: | QA | Karolina Nilsson - KARNIL | 2022-04-12 - 15:48 |
| Approved: | DD | Diana Tieger - DIATIE | 2022-04-14 - 08:06 |
| Released: | QA | Carolina Johansson - SEHRBJNC | 2022-05-19 - 15:07 |

This document has been electronically signed by the persons above.

Provox® LaryButton™



Product description:

Provox LaryButton is delivered single packed, non-sterile, ready for use. The goal is to create a self-retaining, comfortable and airtight fit between the Provox LaryButton and the tracheostoma.

Product Information

| | | | |
|---|--|-----------------|----|
| Document ID: | PF031-01-TechInfo | Edition: | 06 |
| Manufacturer: | Atos Medical AB Kraftgatan 8 SE-242 35 Hörby,Sweden | | |
| Classification: (EU) MDD 93/42/EEC | IIb (2.1 Rule 5) | | |
| Intended Use: | <p>The Provox LaryButton is a self-retaining holder for devices in the Provox HME System intended for vocal and pulmonary rehabilitation after total laryngectomy.</p> <p>For patients with a shrinking tracheostomas it is also used to maintain the tracheostoma for beathing.</p> <p>The Provox LaryButton is intended for single patient use.</p> | | |
| Use specifications: | <p>Intended medical indication Product for rehabilitation for patients breathing through a tracheostoma.</p> <p>Intended patient population Patients of any age. Cognitive ability, by a clinician judged as sufficient. Manual dexterity, by a clinician judged as sufficient. Not intended for patients with mechanical ventilation. Not intended for patients with a low tidal volume.</p> <p>Intended usage Provox LaryButton is a single patient use device prescribed by a clinician.</p> <p>Intended part of the body/type of tissue applied to or interacted with Tracheostoma.</p> <p>Intended user profile The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.</p> <p>Intended conditions of use Environment: Home use (normal daily environment without any or environmental restrictions regarding temperature, moisture etc.). Outpatient clinic use. Hospital use. Frequency of use: Continuous use. Replacement rate: Max usage for 6 months. Replacement is performed by the patient, clinician or caregiver.</p> | | |
| Contraindications: | Provox LaryButton is not intended to be used by patients that are under any form of mechanical ventilation or have damaged tracheostoma tissue. | | |
| CE Mark: | Yes. Devices are CE-marked. | | |
| GMDN code: | 14093 (Tracheostoma button) | | |
| Sterilization: | Non-sterile | | |
| Raw material: | Silicone | | |
| Latex information: | Not manufactured with natural rubber latex | | |
| Biological origin: | The device is not manufactured with materials derived from human or animal source. | | |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C. | | |

Product Information

| | |
|-------------------------------------|--|
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. |
| Hazardous components: | None |
| Expiration date: | 5 years after manufacturing. |
| Packaging: | Provox LaryButton is packed in a plastic bag of polyethylene. The product and instructions for use are packed in a cardboard box. |

Devices under Basic UDI-DI: 7331791-LTU-0-000-0000-38

| REF | Name | UDI-DI |
|------|-------------------------|----------------|
| 7671 | Provox LaryButton 12/8 | 07331791002694 |
| 7672 | Provox LaryButton 14/8 | 07331791002700 |
| 7673 | Provox LaryButton 16/8 | 07331791002717 |
| 7674 | Provox LaryButton 18/8 | 07331791002724 |
| 7685 | Provox LaryButton 12/18 | 07331791002731 |
| 7686 | Provox LaryButton 14/18 | 07331791002748 |
| 7687 | Provox LaryButton 16/18 | 07331791002755 |
| 7688 | Provox LaryButton 18/18 | 07331791002762 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|-------------------------------|---------------------------|
| Provox BasePlate Adaptor | 7331791-HME-A-000-0003-F5 |
| Provox FreeHands HME Cassette | 7331791-HME-0-000-0003-XJ |
| Provox LaryClip | 7331791-LTU-A-000-0001-JT |
| Provox Micron HME | 7331791-HME-0-000-0002-XF |
| Provox ShowerAid | 7331791-ADH-A-000-0000-U8 |
| Provox Swab | 7331791-GEN-A-000-0002-EC |
| Provox TubeBrush | 7331791-GEN-A-000-0001-E9 |
| Provox TubeHolder | 7331791-GEN-A-000-0000-E6 |
| Provox XtraHME | 7331791-HME-0-000-0000-X9 |

Product Information

| Justification: | Function: | Electronic signature justification: | Signed: Date (yyyy-mm-dd) - Time (hh:mm): |
|----------------|-----------|-------------------------------------|---|
| Issued: | QA | Niki Svensson - NIKSVE | 2022-12-19 - 11:36 |
| Reviewed: | QA | Sofia Thomasson - SOFTHO | 2022-12-19 - 11:42 |
| Approved: | QA | Elin Andersson - ELIAND | 2022-12-22 - 15:47 |
| Released: | QA | Niki Svensson - NIKSVE | 2023-03-16 - 13:32 |

This document has been electronically signed by the persons above.



Product description:

The Sizer Kit is a box which contains samples, (Sizers.) of commercially available Provox LaryButtons. The sizes of these Sizers and actual Provox LaryButtons are the same and are indicated on the products themselves and in the bottom of the outer storage box. Each Sizer in the Sizer Kit is stored in an individual removable polypropylene box. This makes it possible for the prescribing specialist to remove the individual storage boxes with the Sizers from the outer storage box individually. This allows for hygienic handling of both the Sizers and the storage boxes. After each sizing session, the Sizer(s) with its individual storage box(es) must be cleaned, disinfected, dried and steam sterilized according to the accompanying Instructions for cleaning and sterilization. The outer storage box must also be cleaned if contaminated. The Sizer LaryButtons and their individual removable storage boxes are thereafter put back at the appropriate position as indicated in the bottom of the outer storage box

Product Information

Document ID: PF032-01-TechInfo **Edition:** 04

Manufacturer: Atos Medical AB
Kraftgatan 8
SE-242 35 Hörby, Sweden

**Classification:
(MDD 93/42/EEC)** IIa (2.1 Rule 5)

Intended Use: The Provox® LaryButton Sizer Kit is intended for use by the prescribing clinician to determine the size(s) of LaryButton that should be prescribed to the patient. The Sizer Kit should be used only by a prescribing clinician who has read the LaryButton Manual. A copy of that manual comes with the Sizer Kit. It can also be viewed on the Internet at www.atosmedical.com. The Sizer LaryButtons are intended for the sizing procedure only. After the correct size(s) have been determined a new LaryButton(s) shall be prescribed to the patient for actual use.

Use specifications: **Intended medical condition**
Laryngectomized patient.

Intended patient population

Gender: Male and female.
Age: Typical average age for a laryngectomy is 65 years.

Intended usage

The Sizer LaryButtons are intended for the sizing procedure only.

Intended part of the body/type of tissue applied to or interacted with
Neck

Intended user profile

Prescribing clinician.

Intended conditions of use

Only to be used in clinical environment.

Contraindications: The Sizer Kit in itself does not have specific contraindications. Do not use the Provox LaryButton, or use it only with special care, in cases of tracheostoma tissue problems such as damaged mucous membrane, granulation tissue formation, and vulnerability with a higher tendency to bleed. The Provox LaryButton may be contraindicated for patients with bleeding disorders or undergoing anticoagulant treatment.

CE Mark: Yes. Device is CE-marked.

GMDN code: 14093 (Tracheostomy button)

Sterilization: Non-sterile, steam sterilizable

Raw material: Silicone, polypropylene

Latex information: Not manufactured with natural rubber latex

Biological origin: The device is not manufactured with materials derived from human or animal source.



Product Information

Handling and storage:

Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.

Waste handling and disposal:

Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.

Hazardous components:

None

Expiration date:

5 years after manufacturing.

Packaging:

Provox LaryButton Sizer Kit is single packed in a tamper-proof plastic bag made of polypropylene together with one IFU for the product, one IFU for Provox LaryButton and one IFU for cleaning and sterilization.

Document No: 10000038474 Edition: 04 Release date: 2023-03-16

Released

Product Information

Devices under Basic UDI-DI: 7331791-LTU-0-000-0001-3B

| REF | Name | UDI-DI |
|------|-----------------------------|---------------|
| 7690 | Provox LaryButton Sizer Kit | 7331791002779 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|-------|--------------|
| N/A | N/A |

| Justification: | Function: | Electronic signature justification: | Signed: Date (yyyy-mm-dd) - Time (hh:mm): |
|----------------|-----------|-------------------------------------|---|
| Issued: | DD | Pontus Eklund - X-PONEKL | 2020-04-21 - 08:24 |
| Reviewed: | DD | Jon Berg - JONBER | 2020-04-21 - 08:58 |
| Approved: | DD | Fredrik Calais - FRECAL | 2020-04-21 - 16:07 |
| Released: | DD | Pontus Eklund - X-PONEKL | 2020-10-28 - 16:39 |

This document has been electronically signed by the persons above.

Product Information

Provox® Fenestration Punch



Product description:

The Fenestration Punch is used for making small fenestrations in a Provox LaryTube at desired locations.

The Fenestration Punch is made of polypropylene, stainless steel and silicone and is used for making small fenestrations in a Provox LaryTube. This is done when the Provox LaryTube is intended to be used in combination with a voice prosthesis.

Document No: 10000038476 Edition: 06 Release date: 2020-10-28

Released

Product Information

| | | | |
|--|--|-----------------|----|
| Document ID: | PF037-01-TechInfo | Edition: | 06 |
| Manufacturer: | Atos Medical AB Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden | | |
| Classification: (EU) 2017/745 | Class I (1.1 Rule 1) | | |
| Intended Use: | The Fenestration Punch is used for making small fenestrations in a Provox LaryTube at desired locations. | | |
| CE Mark: | Yes, the devices are CE marked. | | |
| GMDN code: | 38792 (Basic tracheostomy tube, reusable) | | |
| Sterilization: | Non-sterile | | |
| Raw material: | Stainless Steel, Plastic, Silicone | | |
| Latex information: | Not manufactured with natural rubber latex | | |
| Biological origin: | The device is not manufactured with materials derived from human or animal source. | | |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C. | | |
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. | | |
| Hazardous components: | None | | |
| Expiration date: | 5 years after manufacturing. | | |
| Packaging: | The Fenestration Punch is single-packed in a plastic bag | | |

Product Information

Devices under Basic UDI-DI: 7331791-LTU-A-000-0000-JQ

| REF | Name | UDI-DI |
|------|--------------------------|----------------|
| 7654 | Provox FenestrationPunch | 07331791002632 |

Atos Medical AB Compatible products:

| Range | BASIC UDI-DI |
|-----------------|---------------------------|
| Provox LaryTube | 7331791-LTU-0-000-0002-3E |

| Justification: | Function: | Electronic signature justification: | Signed: Date (yyyy-mm-dd) - Time (hh:mm): |
|----------------|-----------|-------------------------------------|---|
| Issued: | DD | Pontus Eklund - X-PONEKL | 2020-04-20 - 13:35 |
| Reviewed: | DD | Jon Berg - JONBER | 2020-04-20 - 14:57 |
| Approved: | DD | Fredrik Calais - FRECAL | 2020-04-20 - 16:59 |
| Released: | DD | Pontus Eklund - X-PONEKL | 2020-10-28 - 16:35 |

This document has been electronically signed by the persons above.

Product Information

Provox® TubeBrush



Product description:

The Provox TubeBrush is used for cleaning of the Provox LaryTube and Provox LaryButton ex situ. The Provox TubeBrush is packed 6 pieces in a plastic bag. It is available in two different models with outer diameter 8 mm or 12 mm.

Document No: 10000035860 Edition: 09 Release date: 2020-10-28

Released

Product Information

| | | | |
|--|--|-----------------|----|
| Document ID: | PF052-01-TechInfo | Edition: | 09 |
| Manufacturer: | Atos Medical AB Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden | | |
| Classification: (EU) 2017/745 | Class I, Rule 1 | | |
| Intended Use: | The Provox TubeBrush is used for cleaning of the Provox LaryTube and Provox LaryButton ex situ. | | |
| CE Mark: | Yes, the devices are CE marked. | | |
| GMDN code: | 34883 (Airway device, cleaning brush, noninvasive). | | |
| Sterilization: | Non-Sterile | | |
| Raw material: | ABS, Stainless Steel, PBT and Cotton. | | |
| Latex information: | Not manufactured with natural rubber latex | | |
| Biological origin: | The device is not manufactured with materials derived from human or animal source. | | |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C. | | |
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. | | |
| Hazardous components: | None | | |
| Expiration date: | 3 years after manufacturing. | | |
| Packaging: | 6 pieces Provox TubeBrush are packed in a tamperproof plastic bag together with Instructions for Use. | | |

Product Information

Devices under Basic UDI-DI: 7331791-GEN-A-000-0001-E9

| REF | Name | UDI-DI |
|------|------------------------|---------------|
| 7660 | Provox TubeBrush 8 mm | 7331791002656 |
| 7661 | Provox TubeBrush 12 mm | 7331791002663 |

Atos Medical AB Compatible products:

| Range | BASIC UDI-DI |
|-----------------------------|---------------------------|
| Provox LaryTube | 7331791-LTU-0-000-0002-3E |
| Provox LaryButton | 7331791-LTU-0-000-0000-38 |
| Provox LaryButton Sizer Kit | 7331791-LTU-0-000-0001-3B |

| Justification: | Function: | Electronic signature justification: | Signed: Date (yyyy-mm-dd) - Time (hh:mm): |
|----------------|-----------|-------------------------------------|---|
| Issued: | QA | Sara Dahl - X-SARDAH | 2021-11-11 - 18:58 |
| Reviewed: | QA | John Wennborg - JOHWEN | 2021-11-16 - 13:32 |
| Approved: | DD | Diana Tieger - DIATIE | 2021-11-16 - 16:11 |
| Released: | QA | Sara Dahl - X-SARDAH | 2021-12-10 - 09:53 |

This document has been electronically signed by the persons above.

Provox® TubeHolder



Product description:

The Provox TubeHolder has been developed for use with the Provox LaryTube and Provox LaryButton. The integrated clip connectors allow for optimal fit to the wings of the Provox LaryTube and LaryButton, which reduces the physical stress on the soft silicone material.

Product Information

| | | | |
|--------------------------------------|--|-----------------|----|
| Document ID: | PF053-01-TechInfo | Edition: | 07 |
| Manufacturer: | Atos Medical AB Kraftgatan 8 SE-242 35 Hörby, Sweden | | |
| Classification: (EU) 2017/745 | Class I (Rule 1) | | |
| Intended Use: | The Provox TubeHolder is used for extra support for Provox LaryButton and Provox LaryTube. It goes around the neck of the user and the ends are attached to the "ears" of the LaryTube/LaryButton. The Tubeholder is adjustable in length using a Velcro® connection and allows the user to cut the band to suitable length. | | |
| Use specifications: | <p><i>Intended medical indication:</i> Patients breathing through a tracheostoma.</p> <p><i>Intended patient population:</i> Patients of any age. Cognitive ability, by a clinician judged as sufficient. Manual dexterity, by a clinician judged as sufficient.</p> <p><i>Intended usage:</i> Single use. Intended part of the body/type of tissue applied to or interacted with: The device will contact intact skin on the neck.</p> <p><i>Intended user profile:</i> The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.</p> <p><i>Intended conditions of use:</i> Environment: Home use (normal daily environments without any environmental restrictions regarding temperature, moisture etc.). Hospital use. Frequency of use: Continuous use.</p> | | |
| Contraindications: | None. | | |
| CE Mark: | Yes. Devices are CE-marked. | | |
| GMDN code: | 63438 (Tracheostomy tube neck holder, single-use) | | |
| Sterilization: | Non-Sterile | | |
| Raw material: | Tricot textile, Polyurethane (PUR) foam, Polyamide (PA). | | |
| Latex information: | Not manufactured with natural rubber latex | | |
| Biological origin: | The device is not manufactured with materials derived from human or animal source. | | |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C. | | |
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. | | |

Product Information

Hazardous components: None.

Expiration date: 5 years after manufacturing.

Packaging: Single packed together with IFU in a plastic bag.

Devices under Basic UDI-DI: 7331791-GEN-A-000-0000-E6

| REF | Name | UDI-DI |
|------|-------------------|----------------|
| 7668 | Provox TubeHolder | 07331791002670 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|-----------------------------|---------------------------|
| Provox LaryTube | 7331791-LTU-0-000-0002-3E |
| Provox LaryButton | 7331791-LTU-0-000-0000-38 |
| Provox LaryButton Sizer Kit | 7331791-LTU-0-000-0001-3B |
| Provox Life LaryTube | 7331791-LTU-0-000-0004-3L |
| Provox Life LaryButton | 7331791-LTU-0-000-0005-3P |

| Justification: | Function: | Electronic signature justification: | Signed: Date (yyyy-mm-dd) - Time (hh:mm): |
|----------------|-----------|-------------------------------------|---|
| Issued: | QA | Abdallah Almashharawi - ABDALM | 2022-07-25 - 08:36 |
| Reviewed: | QA | Karolina Nilsson - KARNIL | 2022-07-25 - 09:15 |
| Approved: | DD | Peter Sundsten - PETSUN | 2022-07-27 - 08:09 |
| Released: | QA | Abdallah Almashharawi - ABDALM | 2022-07-29 - 09:09 |

This document has been electronically signed by the persons above.

Product Information

Provox® LaryClip™



Product description:

The Provox LaryClip consists of a square adhesive base and a hook-and-loop clip that allows for optimal fit to the wings of the Provox LaryButton and LaryTube. When the adhesive Base is attached to the skin at both sides of the stoma and is eventually removed due to loss of its stickiness, the Clip part can be removed and re-attached as needed.

Product Information

| | | | |
|--------------------------------------|---|-----------------|----|
| Document ID: | PF061-01-TechInfo | Edition: | 08 |
| Manufacturer: | Atos Medical AB Kraftgatan 8 SE-242 35 Hörby,Sweden | | |
| Classification: (EU) 2017/745 | Class I (1.1, Rule 1) | | |
| Intended Use: | The Provox LaryClip is used for extra support for LaryButton and LaryTube. The product consists of two parts, one that is attached to the patients' skin on each side of the stoma and the other part is attached to the LaryButton or the LaryTube. The two parts are then connected by Velcro. | | |
| Use specifications: | <p>Intended medical indication: Product for rehabilitation for patients breathing through a tracheostoma.</p> <p>Intended patient population: Patients of any age. Cognitive ability, by a clinician judged as sufficient. Manual dexterity, by a clinician judged as sufficient.</p> <p>Intended usage: Single use, Over-the-counter</p> <p>Intended part of the body/type of tissue applied to or interacted with: The device will contact intact skin (neck).</p> <p>Intended user profile: The product is supposed to be handled by the patient but is also handled by physicians, trained nurses, SLPs, clinicians and caregivers.</p> <p>Intended conditions of use: Environment: Home use (normal daily environment without any or environmental restrictions regarding temperature, moisture etc.). Outpatient clinic use. Hospital use. Frequency of use: Continuous use. Replacement is performed by the patient, clinician or caregiver.</p> | | |
| Contraindications: | There are no known contraindications. | | |
| CE Mark: | Yes. Devices are CE-marked. | | |
| GMDN code: | 35752 (Tracheostomy tube neck holder, reusable) | | |
| Sterilization: | Non-sterile | | |
| Raw material: | LaryClip Base: Polyethylene (PE), Acrylic Adhesive, velcro LaryClip: Knitted fabric, Polyamide (PA) | | |
| Latex information: | Not manufactured with natural rubber latex. | | |
| Biological origin: | The device is not manufactured with materials derived from human or animal source. | | |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C. | | |

Product Information

| | |
|-------------------------------------|--|
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. |
| Hazardous components: | None. |
| Expiration date: | 3 years after manufacturing. |
| Packaging: | One package consists of 8 pcs of LaryClip and 40 pcs of LaryClip Base. They are packed together with instruction for use in a cardboard box. |

Devices under Basic UDI-DI: 7331791-LTU-A-000-0001-JT

| REF | Name | UDI-DI |
|------|-----------------|----------------|
| 7669 | Provox LaryClip | 07331791002687 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|-----------------------------|---------------------------|
| Provox LaryTube | 7331791-LTU-0-000-0002-3E |
| Provox LaryButton | 7331791-LTU-0-000-0000-38 |
| Provox LaryButton Sizer Kit | 7331791-LTU-0-000-0001-3B |
| Provox Vega Plug 17 | 7331791-VPS-A-000-0004-RU |
| Provox Vega Plug 20 | 7331791-VPS-A-000-0004-RU |
| Provox Vega Plug 22.5 | 7331791-VPS-A-000-0004-RU |

Provox® LaryTube™ Sizer Kit



Product description:

The Sizer Kit is a box which contains samples ("sizers") of a variety of commercially available Provox LaryTubes. The sizes of these Sizers and actual Provox LaryTubes are the same. The size is indicated on the products and both diameter and length are indicated on the chart inside the box. Each sizer in the Sizer Kit is stored in an individual removable polypropylene box. This makes it possible for the prescribing specialist to remove the individual storage boxes with the Sizers from the outer storage box individually. This allows for hygienic handling of both the Sizer(s) and the storage box. After each sizing session, the Sizer(s) with its individual storage box must be cleaned, disinfected, dried and steam sterilized according to the accompanying "instructions for cleaning and sterilization".

| | | | |
|---------------------------------------|---|-----------------|-----|
| Document ID: | PF062-01-TechInfo | Edition: | 2.0 |
| Manufacturer: | Atos Medical AB Kraftgatan 8 SE-242 35 Hörby, Sweden | | |
| Classification: (EU) 93/42/EEC | IIa (Rule 5) | | |
| Intended Use: | The Provox LaryTube Sizer Kit is intended for use by the prescribing specialist to determine the size(s) of LaryTube that should be prescribed to the patient. The Sizer Kit should be used only by a prescribing specialist who has read the LaryTube Manual. A copy of that manual comes with the Sizer Kit. It can also be viewed on the Internet at www.atosmedical.com . The Sizer LaryTubes are intended for the sizing procedure only. After the correct size(s) have been determined, new LaryTube(s) shall be given to the patient for use. | | |
| Use specifications: | Intended medical condition Laryngectomized patient. | | |
| | Intended patient population Gender: Male and female. Age: Typical average age for a laryngectomy is 65 years. | | |
| | Intended usage The Sizer LaryTubes are intended for the sizing procedure only. | | |
| | Intended part of the body/type of tissue applied to or interacted with Neck | | |
| | Intended user profile Prescribing clinician. | | |
| | Intended conditions of use Only to be used in clinical environment. | | |
| Contraindications: | The Sizer Kit in itself does not have specific contraindications. The Provox LaryTubes contained in the LaryTube Sizer Kit are not intended for patients requiring mechanical ventilation. | | |
| CE Mark: | Yes. Device is CE-marked. | | |
| GMDN code: | 12292 (Laryngectomy tube) | | |
| Sterilization: | Non-sterile, steam sterilizable. | | |
| Raw material: | Silicone, Polypropylene. | | |
| Latex information: | Not manufactured with natural rubber latex. | | |
| Biological origin: | The device is not manufactured with materials derived from human or animal source. | | |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2 °C – 42 °C. | | |

Product Information

Waste handling and disposal:

Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.

Hazardous components:

None.

Expiration date:

5 years after manufacturing.

Packaging:

Provox LaryTube Sizer Kit is single packed in a tamper-proof plastic bag together with a manual for the product, instructions for sterilization and a manual for the Provox LaryTube.

Devices under Basic UDI-DI: 7331791-LTU-0-000-0003-3H

| REF | Name | UDI-DI |
|------|---------------------------|----------------|
| 7648 | Provox LaryTube Sizer Kit | 07331791005329 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|-------|--------------|
| N/A | N/A |

Document Approvals
Approved Date: 2023-10-23


| | |
|---|---|
| Task: Approval Task Verdict: Approve | ABDALM Abdallah Almashharawi, Sustaining Engineer (abdallah.almashharawi-atosmedical@coloplast.com) Issuer 16-Oct-2023 07:34:12 GMT+0000 |
|---|---|

| | |
|--|---|
| Task: Final Approval Verdict: Approve | ELIAND Elin Andersson, Associate Design Control & Usability Specialist (elin.andersson-atosmedical@coloplast.com) Technical / Specialist 23-Oct-2023 09:09:38 GMT+0000 |
|--|---|

| Justification: | Function: | Electronic signature justification: | Signed: Date (yyyy-mm-dd) - Time (hh:mm): |
|----------------|-----------|-------------------------------------|---|
| Issued: | QA | Sara Dahl - X-SARDAH | 2021-11-11 - 18:16 |
| Reviewed: | QA | John Wennborg - JOHWEN | 2021-11-16 - 13:34 |
| Approved: | DD | Diana Tieger - DIATIE | 2021-11-16 - 16:12 |
| Released: | QA | Sara Dahl - X-SARDAH | 2021-12-10 - 09:51 |

This document has been electronically signed by the persons above.

Provox® Swab



Product description:

The Provox Swab is a foam attached to a polymer stick handle.

Product Information

| | | | |
|--------------------------------------|--|-----------------|----|
| Document ID: | PF085-01-TechInfo | Edition: | 06 |
| Manufacturer: | Atos Medical AB Kraftgatan 8 SE-242 35 Hörby, Sweden | | |
| Classification: (EU) 2017/745 | Class I, Rule 1 | | |
| Intended Use: | Provox Swab is a single use swab for ex-situ cleaning of Provox LaryTube, Provox LaryButton and tracheostomy inner tubes. | | |
| Use specifications: | <p><i>Intended medical indication:</i> Product for laryngectomized or tracheostomized patients, and/or their caregivers, using Provox LaryTube, Provox LaryButton or double lumen tracheostomy tube, that requires regular cleaning ex-situ.</p> <p><i>Intended patient population:</i> Male and female, laryngectomized or tracheostomized patients.</p> <p><i>Intended usage:</i> Single patient use, swabs should be discarded after use.</p> <p><i>Intended part of the body/type of tissue applied to or interacted with:</i> N/A, cleaning will be performed ex-situ.</p> <p><i>Intended user profile:</i> Patient, clinician, caregiver.</p> <p><i>Intended conditions of use:</i> Normal daily environment without any hygienic or environmental restrictions regarding temperature, moisture etc.</p> | | |
| Contraindications: | No identified or known contraindications. | | |
| CE Mark: | Yes. Devices are CE-marked. | | |
| GMDN code: | 62956 (Airway device cleaning utensil, noninvasive, single-use) | | |
| Sterilization: | Non-Sterile | | |
| Raw material: | Polypropylene (stick handle) and Polyurethane, reticulated foam (foam mitt). | | |
| Latex information: | Not manufactured with natural rubber latex. | | |
| Biological origin: | The device is not manufactured with materials derived from human or animal source. | | |
| Handling and storage: | Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C. | | |
| Waste handling and disposal: | Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard. | | |
| Hazardous components: | None | | |
| Expiration date: | 3 years after manufacturing. | | |

Product Information

Packaging: 50 pcs per package. Devices are packed in plastic bags made of polyethylene and packed together in a cardboard box with printed instructions for use.
Swab Medium is also available as 10pcs, packed in plastic bags with instructions for use printed on the label.

Devices under Basic UDI-DI: 7331791-GEN-A-000-0002-EC

| REF | Name | UDI-DI |
|------|--------------------------|----------------|
| 8250 | Provox Swabs Small | 07331791011412 |
| 8251 | Provox Swab Medium | 07331791011429 |
| 8252 | Provox Swab Large | 07331791011436 |
| 8258 | Provox Swab XtraLarge | 07331791012730 |
| 8083 | Provox Swab Medium 10pcs | 07331791016028 |

Atos Medical AB compatible products:

| Range | BASIC UDI-DI |
|------------------------|---------------------------|
| Provox LaryTube | 7331791-LTU-0-000-0002-3E |
| Provox LaryButton | 7331791-LTU-0-000-0000-38 |
| Provox Life LaryTube | 7331791-LTU-0-000-0004-3L |
| Provox Life LaryButton | 7331791-LTU-0-000-0005-3P |